



ORION INCORPORATION LIMITED

Flat A, 8/F, Block A, 25-31 Kwai Fung Crescent,
Kwai Chung, N.T., Hong Kong

Tel: 852-26159533

Fax: 852-24804234

www.orionzero.hk

GFP Yellow gold electroplating make up

Each unit is for one liter plating solution make up

Prepare: one liter pure distilled water, one liter beaker, heater, glass rod

- 1) Clean beaker, glass rod thoroughly and rinse with pure distilled water
- 2) Fill 800ml pure distilled water into beaker, heat up to 80-85°C, turn off the heat
- 3) Dissolve the small bag of powder completely in the beaker.
- 4) Keep on stir in the big bag of powder while solution temperature around 70°C, until totally dissolved (will take sometime)
- 5) When solution is clear, remove from heat and let it cool down. It is ready to use.
- 6) Solution is ready to use. Please follow operation instructions.



GFP 5N 18K Rose Gold



GFP 2N 18K Yellow Gold



GFP 24K 24K Yellow Gold



GFP 1N 14K Yellow Gold